




Believe in Innovation

20" OF GENIUS



- Results highly correlated to the Westergren reference method
- In accordance with CLSI requirements*
- 1.000 readings for each sample per test
- Use of the same EDTA tube coming from the CBC
- 175 µl of blood
- Self-cleaning without contamination
- Self-checking of reading sensors and physical obstructions
- Data transmission to LIS
- Efficiency and precision
- Standardization and reproducibility
- Certainty of results

EACH SAMPLE IS READ 1000 TIMES IN 20 SECONDS



TEST 1 is the only analyzer capable to give ESR results in 20 seconds by measuring the kinetic of red cells aggregation.

TEST 1 CAPILLARY PHOTOMETRY TECHNOLOGY overcomes the variables and limitations of the sedimentation method also listed in the CLSI document.

SEDIMENTATION ESR	CAPILLARY PHOTOMETRY
Temperature variability	Temperature control 37°C
Hematocrit influence	Independent from Hematocrit value
Dilution problems using Sodium Citrate	No dilution, use of EDTA tube
Inadequate materials and pipettes variability	Use of the same capillary for all samples
Vibrations and pipettes verticality	No influence of vibration or other external factors
Not standardized sample mixing	Automated mixing step
Controls and Calibrators lack	Latex Controls and Calibrators available
Poor reproducibility	High reproducibility

New ESR system generation

SCIENTIFIC PUBLICATIONS

1. CH Cho, GJ Park, YJ Cha, HK Kim, DH Kim, Hoongho, JB Seo, JS Jung, S Jung, HS Cho, DS Lee, H Cho (Department of Laboratory Medicine, University of Seoul College of Medicine and Gangneung Sangji Hospital, Gangneung, Korea). "Erythrocyte Sedimentation Rate Measurements by TEST 1 Better Reflect Inflammation Than Do Those by the Westergren Method in Patients With Malnutrition, Autoimmune Disease, or Infection". *Am J Clin Pathol*. 2008 Feb;119(2):189-94.

2. E. Federo, G. Cipriani, G. Viooli, R. Fontoreno, S. Giommi, I. Costante (Clinical Pathology, Regina Elena Cancer Institute, ICI, Rome - Italy). "Capillary Erythrocyte Sedimentation Rate (ESR) in oncological patients: low hematocrit affects and sample collection optimization in a certified quality system laboratory". *SIBSC National Congress*. 28-31 October 2008, Rimini, Italy.

3. R. Pipipi, E. Piva, B. Ribonari, F. Tosato, M. Pibiani, (Dept. of Laboratory Medicine, Pavia University School of Medicine, Pavia, Italy) "The Erythrocyte Sedimentation Rate (ESR): an old test with new contents". *SIBSC National Congress*. 28-31 October 2008, Rimini, Italy.

4. J. Rile, J. Diamentino, N. Carini, F. Malici (Clinical Pathology Department, IEO Coimbra, Francisco Gentil, EPE, Portugal). Erythrocyte sedimentation rate in blood: a comparison of the Test 1 ESR system with the CGH reference method. *Clinical Chemistry and Laboratory Medicine*. 2007 June; 45, Special Supplement, p. S116, M277.

5. E. Piva, R. Pipipi, V. Tomponi, M. Pibiani (Dept. of Laboratory Medicine, University Hospital, Padova, Italy) "A new turboanalytic standard to improve the quality assurance of the erythrocyte sedimentation rate measurement". *Clinical Biochemistry* 2007 Apr; 42(7):491-5. Epub 2007 Jan 8.

6. Arslan, N. Arslan, N. (Biochemistry Department, Basirent University, Ankara, Turkey) "Comparison of the erythrocyte sedimentation rate measured by the Micro Test 1 sedimentation analyzer and the conventional Westergren method". *Ann Saudi Med* 2007; 27(6): 382-385.

7. Li J, Yi WB, Chen G, Fang SF (Shen Dept of Clinical Laboratory, Sixth People's Hospital Affiliated to Shanghai Jiaotong University, Shanghai 20023, China). "Evaluation of the Microtest 1 ESR analyzer and investigation of the reference value". *Chin J Lab Med* 2007; 30(2): 103-104 (article in Chinese).

8. S. Ockers, HS Abbas, L. Donnee, M. Quibben (Clinical Biochemistry Unit, Medical Faculty, Central Laboratory, Akdeniz University, Antalya, Turkey) "Comparison of TEST 1 with SRS 100 and CGH reference method for the measurement of the length of sedimentation reaction in blood". *Clin Chem Lab Med* 2008; 46(4):607-12.

9. NE Ajaja, AJ Bakker, GA van den Berg (Slingsted KCL, Dept. of Clinical Chemistry, Leeuwarden, The Netherlands). "Determination of the length of sedimentation reaction in blood using the Test 1 system, comparison with the Sedimetric 100 method, sedimentation boronog levels and the influence of M protein". *Clin Chem Lab Med* 2008, 44(7): 904-909.

10. Y. Kawase, N. Inada, S. Ito, S. Makino, M. Miyake (Clinical Test Section of Eiji Hospital, Fingag Line Co. and Dept. of Clinical Pathology of Juntendo University, Evaluation for ESR automated measuring instrument with EDTA³, 36th Japan Society for Clinical Laboratory Automation, 30 September 2008, Japan).

11. B. Rissak, P. Chai, C. Maia, J. Jakubina (Servicio Laboratorio Clínico, Hospital Universitario de Cangas, Universidad Comarcal de Santiago de Equador que realiza VHS, Test 1 y Venenec³, VII Congreso Chileno de Tecnología Médica, 20-22 October 2008, Santiago, Chile (article in Spanish)).

12. M. Pibiani, P. D'Abate, V. Tomponi, E. Piva, M. Buvarello, M. Sarconi, (Dept. of Laboratory Medicine, University Hospital, Padova, Italy). "Validazione della tecnica in ed individualità della velocità d'Entrosedimentazione". 30th SIBSC, 8-11 June 2008, Padova, Italy (article in Italian).

13. E. Melik, M. Pipipi, P. Lorenz (Klinički center Ljubljana, Klinični inštitut za klinično kemijo in biokemijo "Novi nudi merjenje in študije sedimentacije eritrocitov v endokrinološki testni kitovi". 2. kongres hematologije in transfuziologije Slovenije in Medicinarsko Ljubljano, 23-24 April 2008, Portoroz, Slovenia (article in Slovenian)).

14. H. Chiva Alonso, M. Girvent Moreno, MT Robles Bello, V. Balcells Arco, O. Vidal JM (Laboratorio Hospital San Vicente y Area Sanitaria, H. Algora, Sagunto, Camino De Melisera Para La Delimitación De VEG, G. Reparaciones Sobre La Fase Preanalítica", Generalitat Valenciana - Conselleria De Sanitat for Valencia Government - MOH), Spain 2008 (paper in Spanish).

15. P. Galloni, "Quality and Automation in the Determination of the Erythrocyte Sedimentation Rate", Symposium 248, 2004 World Congress of Pathology & Laboratory Medicine, 30 August - September 2004, Busan, Korea.

16. M. Nicolò, E. Lanza, A. Maresco, (Laboratory of Clinical Chemistry and Hematological Analysis, Ospedale Civile Maggiore, Verona, Italy) "Integrated Haematology and Coagulation Laboratory". Poster, EuroMedLab Congress, 1-5 June 2003, Barcelona, Spain.

17. M. Pibiani (Dept. of Laboratory Medicine, University Hospital, Padova, Italy). "Erythrocyte Sedimentation Rate: Innovative Techniques for an Obsolete Test?". *Clinical Chemistry and Laboratory Medicine*, 2003, 41 (2): 115-116.

18. A. Barone, M. Mancini, G. Barone (Dept. of Haematology, HUGO, "Spagn de la Vico", Azienda S. SIMEASA, School of Medicine, University of Milano, Spini), "Determination of the Length of Sedimentation Reaction in Blood: a Comparison of the Test 1 ESR System with the CGH Reference Method and the Sedimetric³". *Clinical Chemistry and Laboratory Medicine* 2003, 41 (2).

19. D. Calzavara, S. Capozzo, F. Caciulli, M. Cariti, G. Soffici (Clin. Chem. e Hematol. Lab., San Bonifacio Hospital, Verona, Italy). "Internal Quality Control for Erythrocyte Sedimentation Rate Measured Test 1 Analyzer". *Clinical Laboratory* 2002, 49:459-462.

20. E. Nicosi, (Istituto Nazionale di Tecnologia, Interlab). "Comparison of the Westergren method versus the TEST1 technique for determining the Erythrocyte Sedimentation Rate". May 2002; private communication.

21. BH Lee, J. Choi, MS Gwe, KK Lee, H Park (Dept. of Laboratory Medicine, Kyungju Samsung Hospital, Sungju-gun, Kyungpook University School of Medicine, Seoul, Korea). "Stem Evaluation and Reference Range Assessment of TEST1 for the Automated Erythrocyte Sedimentation Rate". *Journal of Clinical Pathology and Quality Control*, Vol. 24, No. 1, 2002 (article in Korean).

22. E. Piva, P. Passari, M. Pibiani (Dept. of Laboratory Medicine, University Hospital, Padova, Italy). "Determination of the length of sedimentation reaction (erythrocyte sedimentation rate) in non-anticoagulated blood with the Microtest 1". *Clin Chem Lab Med* 2003; 41:677-737.

23. M. Pibiani, E. Piva (Dept. of Laboratory Medicine, University Hospital, Padova, Italy). "Erythrocyte Sedimentation Rate: Use of Fresh Blood for Quality Control". *American Journal of Clinical Pathology*, 2002, 111:621-626.

24. D. Smith D. Spedding, (Clyde Balming Diagnostics, New Zealand). "Evaluation of Agreement between the TEST1 and Stansed Automated ESR Analyzers". November 2001; private communication.

25. D. Calzavara, S. Capozzo, M. Cariti, F. Caciulli, G. Soffici (Clin. Chem. e Hematol. Lab., San Bonifacio Hospital, Verona, Italy). "Internal Quality Control for Erythrocyte Sedimentation Rate (ESR) measured by TEST 1 Analyzer". *Clinical Chemistry*, June 2001, 47: 162.

26. E. Piva, M. Sarconi, G. Barone, M. Pibiani (Dept. of Laboratory Medicine, University Hospital, Padova, Italy). "Length of Sedimentation Reaction in Unclotted Blood (Erythrocyte Sedimentation Rate): Variation with Sex and Age and Reference Limits". *Clinical Chemistry and Laboratory Medicine*, May 2001, 39: 451-454.

27. D. de Jongh, J. Siewkwaning, J. Sings, JM Pijls (Dept. Clinical Chemistry, University Hospital, The Netherlands). "Erythrocyte Sedimentation Rate by Test 1 Analyzer". *Clinical Chemistry*, June 2000, 46: 881-882.

28. M. Pibiani, S. De Tosi, M. Sarconi, G. Barone, E. Bolognini (Department of Laboratory Medicine, University Hospital of Padova, Italy). "The TEST 1 automated system: a new method for measuring the erythrocyte sedimentation rate". *Am J Clin Pathol*. 1998 Sep; 110(3):234-40.

29. G. Soffici (Clinical Chemistry and Hematology Laboratory, San Bonifacio Hospital, Verona, Italy). "Nuovo Metodo per la Determinazione della Velocità di Entrosedimentazione (VES)". *Biochimica Clinica*, Vol. 22, N. 6, 1998, p. 300.

30. N. Cava, Z. Abu-Arg, N. Giacconi, F. Biondini, G. Pizzardi, M. Troceni (Dept. of Laboratory Medicine, G. Sassi Hospital, Ancona, Italy). "TEST 1: un Nuovo Metodo per la Determinazione della VES". *Biochimica Clinica*, Vol. 22, N. 6, 1998, p. 300.



ESR in 20 Seconds



* Clinical Laboratory and Standard Institute (formerly NCLSI). "Reference and Selected Procedure for ESR Test; Approved standard- 4th Edition" Vol.20 No.27; December 2004

Rev. 22/03/2011

Believe in Innovation

FULLY AUTOMATED ANALYZER FOR THE DETERMINATION OF THE ERYTHROCYTE SEDIMENTATION RATE



- Results in 20 seconds related to red cells aggregation
- First result available after 5 minutes from analysis start
- No reagents required
- Results expressed in mm/h
- High correlation with the Westergren method
- No influence of low hematocrit levels
- Use of the same CBC tubes
- 175 µl EDTA blood sample per test
- Only 800 µl sample requested in the tube
- Capacity up to 60 samples
- Throughput up to 180 samples/hour
- Latex Calibration & Controls
- Smart cards
- Thermostated at 37°C
- Mixing cycles in accordance with CLSI requirements*
- Internal bar code reader
- Connection to LIS
- Simplified needle replacement
- Thermal printer
- Complete integration in the hematology laboratory



Test 1 THL

Code SI 195.210/THL

Up to 60 samples per session with standard Test 1 racks. Call Blood Counter racks adaptors available.



CE MARKED

FULLY AUTOMATED ANALYZER FOR THE DETERMINATION OF THE ESR WITH EXTERNAL NEEDLE FOR PEDIATRIC SAMPLES AND UNCAPPED TUBES



- Results in 20 seconds related to red cells aggregation
- First result available after 5 minutes from analysis start
- No reagents required
- Results expressed in mm/h
- High correlation with the Westergren method
- No influence of low hematocrit levels
- Use of the same CBC tubes
- 100 µl EDTA blood sample with manual withdrawal and 175 µl with automated withdrawal per test
- Only 100 µl sample requested in the tube for manual withdrawal and 800 µl for automated withdrawal
- Latex Calibration & Controls
- Smart cards
- Thermostated at 37°C
- Mixing cycles in accordance with CLSI requirements*
- LCD touch screen
- User-friendly software
- Connection to LIS
- Simplified needle replacement
- Thermal printer



Code SI R10PN, SI R20PN

CE MARKED



Roller 10PN/20PN

Technical Features

Power supply: 115-230 VAC, 50/60 Hz Power consumption: 100VA Operative temperature: from +10 to +30°C Size: 240x300x400 mm Weight: 14/10 Kg Initial specifications: two RS232 serial ports Optional accessories: external CCD bar code reader

LATEX CONTROLS



THREE LEVELS TO GUARANTEE:

Precision
Accuracy
Repeatability



Latex Control Kit 6 tests
Code SI 305.100-A (Greiner tubes), SI 305.102-A (Sarstedt tubes)



Latex Control Kit 30 tests
Code SI 305.300-A (Greiner tubes), SI 305.302-A (Sarstedt tubes)

Alifax ESR External Quality Programs successfully started in many Countries!

External Quality Evaluation Kit

Specific for Test 1 family instruments
Code SI 305.500-A (Greiner tubes), SI 305.502-A (Sarstedt tubes)

NEW DIRECT LOADING FROM THE ORIGINAL CELL BLOOD COUNTER RACK



Test 1 BCL Code SI 195.220/BCL
Test 1 SDL Code SI 195.230/SDL
Test 1 YDL Code SI 195.240/YDL
Test 1 MDL Code SI 195.250/MDL
Test 1 XDL Code SI 195.260/XDL

Up to 60 samples per session with ALIFAX green plastic racks. Up to 48 samples per session with Beckman Coulter CBC racks.

Up to 40 samples per session with ALIFAX yellow plastic racks. Up to 40 samples per session with Sysmex CBC racks.

Up to 40 samples per session with ALIFAX blue plastic racks. Up to 40 samples per session with Siemens CBC racks.

Up to 40 samples per session with Beckman Coulter CBC racks.

Up to 40 samples per session with Beckman Coulter CBC racks.

Test 1 - THL-BCL-SDL-YDL-MDL-XDL Technical Features
Power supply: 115-230 VAC 10% 50/60 Hz Power consumption: 150 VA max Operative temperature: from +10 to +31 °C Size: 510x600x600 mm Weight: 45 Kg Initial specifications: two RS232 serial ports LCD bar reader: Internal



- Results in 20 seconds related to red cells aggregation
- First result available after 5 minutes from analysis start
- No reagents required
- Results expressed in mm/h
- High correlation with the Westergren method
- No influence of low hematocrit levels
- Use of the same CBC tubes
- Only 800 µl sample requested in the tube
- Latex Calibration & Controls
- Smart cards
- Thermostated at 37°C
- Mixing cycles in accordance with CLSI requirements*
- Connection to LIS
- Simplified needle replacement
- Thermal printer
- Automated washing



Code SI R20-LC

CE MARKED



Roller 20LC

Technical Features

Power supply: 115-230 VAC, 50/60 Hz Power consumption: 50 VA max Operative temperature: from +10 to +30 °C Size: 300x600x600 mm Weight: 23,2 Kg Initial specifications: two RS232 serial ports Optional accessories: external CCD bar code reader

NEW TWO PARAMETERS CARDS FOR TEST 1 FAMILY INSTRUMENTS

Smart CARDS



1000 tests



4000 tests



1000 tests



2000 tests

Environmentally friendly cards save storage and transport costs. Only pay for a test when needed.

