

Samsung Electronics Co., Ltd. inspires the world and shapes the future with transformative ideas and technologies, redefining the worlds of TVs, smartphones, wearable devices, tablets, cameras, digital appliances, printers, medical equipment, network systems, and semiconductor and LED solutions. We are also leading in the Internet of Things space through, among others, our Smart Home and Digital Health initiatives.

We employ 307,000 people across 84 countries with a annual sales of US \$196 billion. To discover more, please visit our official website at [www.samsung.com](http://www.samsung.com) and our official blog at [global.samsungtomorrow.com](http://global.samsungtomorrow.com)

CT-LABGEO IB10 4.5-150805-EN

Thermo Fisher Scientific is the global collaboration partner of Samsung for commercialization and co-development of innovative solutions in the POC testing market. Together we provide excellent expertise in clinical diagnostics to fulfill customer needs and improve patient management.

Samsung Smart Blood Analyzer

IMMUNOASSAY ANALYZER



Samsung LABGEO HC10

HEMATOLOGY ANALYZER



Samsung LABGEO PT10

CLINICAL CHEMISTRY ANALYZER

Samsung  
**LABGEO** IB10

Take Your Laboratory Wherever You Want



**SAMSUNG**

© 2015 Samsung Electronics All rights reserved.  
Samsung Electronics reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.

**SAMSUNG**



# SAVE TIME, SAVE LIVES

## SMART SOLUTION PROVIDER

### Samsung LABGEO IB10

High Performance Portable  
Immunoassay Analyzer

WWW.INSIGNIABIOMEDICAL.COM

### Fast Easy Accurate

The Samsung LABGEO IB10 is a portable, reliable and easy-to-use immunoassay analyzer for *in-vitro* determinations and provides medical professionals with rapid and accurate test results in emergency care settings where near patient testing is critical.



# FAST

# ACCURATE

WWW.INSIGNIABIOMEDICAL.COM

## Fast & Fully Automated

Rapid test results in 20 minutes  
Fully automated centrifugation inside the analyzer

## Whole Blood Use

No sample pretreatment required for plasma separation

## Multiple Analytes Tested in a Single Run

No additional time required to test multiple analytes  
Up to three analytes tested in a single run – One run, Three test results in 20 minutes

# EASY

## Simple & Easy-to-Use

Simple, LCD touch screen interface with clear, step-by-step instruction  
Multiple language support

## Convenient & User Friendly

One step convenient testing method  
Test Discs stable at room temperature for 30-days  
Built-in thermal printer  
Multiple external devices available

## Portable

Easy to move (Compact size, Lightweight)  
Optional battery operation  
Accurate testing when and where it's needed



## Excellent Precision and Correlation

Excellent assay performance assists medical professionals in the diagnosis and treatment of critically ill patients

	Troponin I	CK-MB	Myoglobin	NT-proBNP	D-Dimer	Beta-hCG	TSH
Measurable Range	0.05 – 30 ng/mL	2.0 – 60 ng/mL	30.0 – 500 ng/mL	30 – 5000 pg/mL	100 – 4000 ng/mL (FEU)	4–400 mIU/mL	0.11–120.00 mIU/L
Upper Reference Limit (URL)	0.10 ng/mL (99 %ile)	8.58 ng/mL	99.84 ng/mL	< 75yrs = 125 pg/mL ≥ 75yrs = 450 pg/mL	446.8 ng/mL (FEU)	5.42 mIU/mL	0.25–3.85 mIU/L (euthyroid reference interval)
Precision	Within-run CV	9.4 – 12.6%	9.0 – 11.2%	8.5 – 10.2%	6.7 – 8.6%	12.4–13.4%	5.8–11.2%
	Total CV	10.1 – 13.0 %	10.3 – 11.5%	12.1 – 12.7%	9.0 – 10.9%	6.7 – 9.0%	6.0–11.4%

\* Each laboratory should establish a reference range that represents the patient population that is to be evaluated.

## Method Comparisons Show High Correlation with Central Laboratory Systems

Correlation Coefficients above 0.90 and method comparison slopes near unity compared to central laboratory systems

Figure 1. Method Comparison – LABGEO<sup>®</sup> NT-proBNP and Ortho VISTA<sup>®</sup> NT-proBNP





Figure 2. Method Comparison – LABGEO<sup>®</sup> CHF NT-proBNP and Ortho VISTA<sup>®</sup>/Roche NT-proBNP



The Samsung LABGEO<sup>®10</sup> Analyzer produced rapid and accurate measurement of both Troponin I and NT-proBNP in whole blood and plasma and showed high correlation with the central laboratory system, the Siemens Dimension Vista<sup>™</sup> platform

Table 3. Method Comparison - LABGEO<sup>®</sup> TnI and Siemens Dimension Vista<sup>™</sup> cTnI\*

Comparison	N	Intercept (ng/mL)	Slope	Spearman rank correlation
LABGEO <sup>®</sup> WB vs. VISTA plasma	36	-0.070	0.868	0.911
LABGEO <sup>®</sup> plasma vs. VISTA plasma	36	-0.084	0.787	0.915
LABGEO <sup>®</sup> WB vs. LABGEO <sup>®</sup> plasma	36	0.006	0.993	0.964

Table 4. Method Comparison - LABGEO<sup>®</sup> NT-proBNP and Siemens Dimension Vista<sup>™</sup> NT-proBNP\*

Comparison	N	Intercept (pg/mL)	Slope	Spearman rank correlation
LABGEO <sup>®</sup> WB vs. VISTA plasma	47	-6.6	1.028	0.938
LABGEO <sup>®</sup> plasma vs. VISTA plasma	47	19.1	0.881	0.926
LABGEO <sup>®</sup> WB vs. LABGEO <sup>®</sup> plasma	47	-9.5	1.062	0.978

\*A study was performed by Dr. Steven Kazmierczak at Oregon Health and Science University, OHSU (Portland, Oregon) AACCC 2013 poster Kazmierczak, S. (suppl) pg S2-S310 - B221

## AVAILABLE TEST MENU

Product Image	Product Name	Model #
	Samsung LABGEO® Troponin I Test Troponin I	IVR-IB05
	Samsung LABGEO® CHF Test NT-proBNP	IVR-IB06
	Samsung LABGEO® CHF Plus Test Troponin I / NT-proBNP	IVR-IB07
	Samsung LABGEO® Cardiac 3-in-1 Panel Test Troponin I / CK-MB / Myoglobin	IVR-IB08
	Samsung LABGEO® D-Dimer Test D-Dimer	IVR-IB09
	Samsung LABGEO® SOB Cardiac 3 Panel Test Troponin I / NT-proBNP / D-Dimer	IVR-IB11
	Samsung LABGEO® beta-hCG Test beta-hCG	IVR-IB12
	Samsung LABGEO® TSH Test TSH	IVR-IB15

## Rapid Point-of-Care (POC) Immunoassay



**Troponin I Test**  
Troponin I

For the Quantitative Determination of Troponin I (cTnI) in Lithium Heparin Human Whole Blood and Plasma



**CHF Test**  
NT-proBNP

For the Quantitative Determination of N-terminal pro-brain natriuretic peptide (NT-proBNP) in EDTA and Lithium Heparin Human Whole Blood and Plasma



**CHF Plus Test**  
Troponin I / NT-proBNP

For the Quantitative Determination of Troponin I (cTnI) and N-terminal pro-brain natriuretic peptide (NT-proBNP) in Lithium Heparin Human Whole Blood and Plasma



**Cardiac 3-in-1 Panel Test**  
Troponin I / CK-MB / Myoglobin

For the Quantitative Determination of Troponin I (cTnI), CK-MB, Myoglobin in Lithium Heparin Human Whole Blood and Plasma



**D-Dimer Test**  
D-Dimer

For the Quantitative Determination of D-Dimer in Lithium Heparin and Citrated Human Whole Blood and Plasma



**SOB Cardiac 3 Panel Test**  
Troponin I / NT-proBNP / D-Dimer

For the Quantitative Determination of Troponin I (cTnI), NT-proBNP and D-Dimer in Lithium Heparin Human Whole Blood and Plasma

## SIMPLE ONE STEP PROCEDURE




## KEY FEATURES

<b>Analyzer (Size / Weight)</b>	177(W) x 330(D) x 177(H) mm / 2.4 kg
<b>Cardiac Markers</b>	Troponin I, NT-proBNP, D-Dimer, CK-MB, Myoglobin
<b>Product Menu</b>	Troponin I Test, CHF Test, CHF Plus Test, 3-in-1 Test, D-dimer Test, SOB Cardiac3 Panel Test, Beta-hCG Test, TSH Test
<b>Sample Volume</b>	500µl Whole Blood or Plasma
<b>Options</b>	Battery, Barcode Scanner, WIFI USB Dongle, EQC Disc
<b>IT Solution</b>	Samsung LABGEO Data Manager
<b>Proactive Service</b>	RMS (Remote Management System)

## SMART HEALTHCARE

### Proactive Customer Service

RMS (Remote Management System) provides real-time service through equipment status checks and immediate response to errors from a service team



### IT Solution

Connects to HIS (Hospital Information System), PC, and Mobile phone / device for efficient database management



### Support Multiple Devices

Multiple external devices available (WIFI USB Dongle, USB, Barcode Scanner, and LAN Connectivity)

